

INDEX

A

- Adaptive immunity.....3–4, 6, 8, 9, 186, 216, 232, 244, 373, 422, 423
- Adduct..... 52, 53, 125
- Adhesion 4, 75, 76, 133, 145, 216, 217, 269–271, 275–276, 280, 281, 283, 286, 290, 293, 294, 318, 320, 323–324, 350, 452
- Affinity/binding affinity..... 6, 21, 39, 43, 44, 48, 55, 59, 67, 68, 71, 78–84, 86–89, 92, 94, 96, 101, 102, 105, 106, 108, 109, 111, 115, 128, 153, 174, 195, 212, 281, 360, 361, 363, 403, 411–412, 415, 446
- Alkylate 53, 55, 57, 58, 193, 198
- Apoptosis..... 4, 6–8, 11, 14, 16, 17, 77, 153, 159, 160, 185, 186, 193–195, 232, 246, 251, 266, 328, 355, 356, 408, 451–461
- Asialofetuin (ASF)..... 78, 80, 84, 86

B

- Bacteria..... 15, 39–42, 47, 48, 117, 119, 123, 124, 154, 157, 175, 202, 204, 205, 208–210, 212, 213, 226, 232, 234, 237, 244, 329, 335, 372, 382, 422, 424–429, 472
- β -galactose..... 115, 116
- β -galactosidase..... 95, 97, 100, 413
- Binding..... 3, 38, 52, 63–73, 76, 91–103, 105–113, 115–128, 147, 154, 170, 194, 202, 216, 245, 250, 271, 294, 318, 328, 351, 355, 372, 398, 422, 448, 452, 472
- Binding assay..... 100, 101, 245, 385, 417
- Biopolymer 79
- Biotinylation..... 95, 97, 121–122, 125, 271, 280
- Blood group..... 21, 79, 232, 421–423, 427, 429

C

- Carbohydrate-recognition domain (CRD)..... 18–20, 38, 39, 219, 322, 328, 343, 409
- cDNA..... 38, 41, 48, 77, 154–156, 331, 361, 453
- Cell death..... 4, 8, 10, 11, 76, 155, 185–187, 193, 194, 265, 459
- Cell growth..... 16, 75, 76, 148, 451–461
- Cell surface glycans 2, 7, 21, 92–95, 126, 127, 226, 250, 452
- Centrifuge 42, 69, 117, 119, 120, 123, 138–140, 147, 149, 160, 172, 174–176, 180, 189, 190, 192, 202, 203, 205–210, 238, 239, 242, 243, 256–258, 260, 261, 273–275, 279, 289, 298, 301, 335, 376,

- 377, 379–382, 393, 403, 404, 407–409, 434, 440, 442–444, 453, 454, 457–460, 466, 470, 480
- Chromophore..... 78
- Clustered regularly interspaced short palindromic repeats(CRISPR-Cas) 329, 330, 338
- Coding sequence 264, 328, 473, 476, 477
- Column 40, 54, 64, 87, 95, 118, 155, 171, 188, 233, 360, 375, 403, 424, 446
- Comparative glycomics..... 63
- Concanavalin A (ConA)..... 80, 81, 84
- Covalent 52, 67, 78, 83, 94, 105, 111, 153, 182, 357, 414, 445
- CRISPR-Cas. *See* Clustered regularly interspaced short palindromic repeats(CRISPR-Cas)
- Cross-link..... 20, 38, 47, 76, 77, 83, 85, 86, 115, 281, 351, 398, 405, 445, 446, 455
- Cysteine..... 52, 53, 94, 95, 98, 357, 412, 415

D

- Dendritic cell (DC)..... 9, 10, 16, 215–227, 240, 250, 254–255, 260–262
- Digest..... 55, 58, 78, 88, 99, 101, 201, 238, 242, 245, 247, 334, 335, 338, 441, 472, 473
- Dimer..... 17, 18, 20, 54, 55, 95, 158, 159
- Dimerization 17, 18, 20, 52, 54, 55, 415
- Disaccharide 19, 322
- Dissociation constant 63, 67–70, 84, 86, 88, 101
- Dissociation rate..... 80
- Disulfide bond..... 52, 53
- Divalent..... 281
- Domain 18–20, 38, 39, 47, 75–79, 93, 115, 170, 175, 322, 329, 357, 409, 436, 443–444, 452

E

- Elution 40, 42–45, 48, 55, 57, 59, 64, 67, 68, 70, 72, 118, 120–122, 155, 157, 158, 165, 171, 174, 175, 415
- Embryo..... 8, 13, 305, 306, 308–310, 314, 327, 328, 330–338, 399, 405
- Endo- β -galactosidase 95, 97, 100, 413
- Endoplasmic reticulum (ER)..... 51, 52, 60
- Enzymatic 78, 92, 97, 102, 246, 247
- Epithelia(l) 4, 7, 12, 13, 16, 37, 49, 76, 232, 234, 237–241, 245–247, 317–320, 322–324, 344, 366, 374
- Epitope..... 75, 77, 80, 81, 84, 86, 231, 447, 479
- Epsilon binding protein (eBP) 11, 12

Equilibrium dialysis.....	78
ER. <i>See</i> Endoplasmic reticulum (ER)	
Erythrocyte.....	87, 274, 275, 378, 386, 387
Eukaryotic.....	453
Expression.....	5, 37–49, 76, 119, 148, 154, 170, 185, 202, 216, 243, 250, 270, 289, 294, 310, 318, 327–339, 343, 368
F	
FAC. <i>See</i> Frontal affinity chromatography (FAC)	
FITC. <i>See</i> Fluorescein isothiocyanate (FITC)	
Flow cytometric analysis.....	91–103, 122, 143, 194–196, 204–205, 209–210, 245, 257, 350
Fluorescein isothiocyanate (FITC).....	122, 147, 156, 162, 172, 173, 188, 193, 198, 204–206, 209–210, 212, 213, 217, 233, 235, 253, 254, 259, 264, 271–273, 275, 276, 278, 279, 281, 458, 459, 467, 471
Fluorescence labeling.....	64, 95, 96, 98, 99, 203, 205–207, 212
Frontal affinity chromatography (FAC).....	19, 63–73, 106, 116
G	
Galactose.....	20, 65, 375, 379, 387, 389
Galectin	
canine galectin.....	37–39
domain	
C-terminal.....	20, 436, 443
N-terminal.....	20, 436, 443, 452
GBP. <i>See</i> Glycan binding protein (GBP)	
Gene expression.....	310, 328, 329, 356, 453–454, 456–457
Glutathione S-transferase (GST)-galectin	
galectin	
chimera-type.....	5, 328
galectin-1.....	3, 38, 51–61, 71, 76, 94, 106, 125, 153, 169, 193, 292, 294, 306, 343, 371, 397, 429, 432, 451, 465
galectin-2.....	3, 6, 11, 17, 19, 20, 38, 108, 109, 123, 125, 128, 232, 350, 452
galectin-3.....	3, 37, 70, 76, 108, 125, 154, 202, 251, 294, 318, 343, 355–368, 371, 397, 429, 431, 451, 465–482
galectin-4.....	3, 6, 11, 15, 20, 21, 37, 38, 232, 344, 422, 423, 429
galectin-7.....	3, 13, 14, 16, 17, 38, 39, 125, 318, 344, 452, 465
galectin-8.....	3, 6–8, 11, 13, 15, 20, 21, 38, 39, 93, 127, 154, 164, 166, 255, 270, 275, 276, 278, 279, 294, 422, 423, 426, 429
galectin-9.....	3, 6, 8–11, 13, 17, 20–22, 38, 39, 43, 46, 47, 202, 232, 251, 255, 262, 371
galectin-10.....	3, 6, 11, 38, 375
human.....	71, 185–199, 220, 322
PpGalec.....	22
proto-type.....	328
recombinant.....	54, 57, 59, 65, 95, 119, 121–124, 141, 146, 149, 154, 166, 188, 190–192, 223, 226, 227, 271–273, 286–289, 319, 324, 375, 378–379, 386–389, 401–403, 412, 424, 426, 428
sandfly <i>Phlebotomus papatasi</i>	22
tandem-repeat-type.....	3, 18–20, 39, 47, 328
zebrafish.....	13
Glycan binding protein (GBP).....	2, 3, 17, 91, 94, 115, 116, 125, 126, 250, 294
Glycoconjugates.....	92, 94
Glycoprotein.....	2, 8, 38, 75, 76, 78–80, 84, 91, 126, 156, 216, 219, 221, 279, 281–283, 398, 452
Glycoprotein receptor.....	3, 75, 76
Gradient affinity chromatography.....	78
<i>Griffonia simplicifolia</i> lectin II (GS-II).....	96, 99, 100
H	
Hemagglutination.....	80, 87, 378–379, 386–387, 409, 410, 415
inhibition.....	78, 81, 82, 87, 116
High-performance liquid chromatography	
(HPLC).....	55–56, 58–60, 64, 67
Homology.....	38, 328
Human immunodeficiency virus (HIV).....	8, 21, 22, 356, 400–401, 405–408, 413
I	
IFN- γ . <i>See</i> Interferon gamma (IFN- γ)	
Immobilize.....	20, 63, 65–70, 72, 93, 94, 96, 97, 99–102, 105, 112, 270, 271, 275–276, 281, 351, 375, 438, 446
Incubator.....	119, 135, 138, 139, 160, 161, 173, 203, 205, 206, 208, 210, 222, 223, 252, 253, 256, 262, 263, 286–289, 294–296, 298, 299, 307, 309, 312, 313, 321, 324, 333, 334, 336, 359, 376, 377, 380, 381, 383, 384, 399, 401–403, 424, 425, 427, 453, 455–457, 471, 475
Inflammation.....	9, 11, 17, 75, 76, 133, 186, 215–217, 219–220, 231–247, 250, 283, 317
Inhibitory.....	8, 80, 250, 290
Innate immunity.....	2, 8, 21–22, 51, 231, 232, 328, 421–423, 429
Insoluble.....	157
Insoluble precipitate.....	89
Interferon gamma (IFN- γ).....	4, 6, 9, 10, 232, 254, 260, 454, 457
Interleukin-5 (IL-5).....	6, 7, 254, 260
Interleukin-6 (IL-6).....	6
Interleukin-10 (IL-10).....	4–6, 9, 10, 245, 249, 250, 254, 260
Iodoacetamide.....	51–61, 193, 198, 401, 409, 410, 416
Isothermal titration calorimetry (ITC).....	75–89

L

- Lactose 39, 40, 42–48, 54, 55, 57,
59, 95, 96, 98, 100–102, 106, 108, 117, 118, 121,
122, 150, 155, 157–162, 171, 188, 190, 191,
193–197, 296, 299, 319, 321, 322, 324, 379, 387,
388, 390, 401–403, 412, 415, 426, 428, 429, 437,
438, 445, 446, 452, 481
- LEA. *See Lycopersicon esculentum* (tomato)
agglutinin (LEA)
- Lectin 68, 73, 75–89, 99, 106, 202,
231, 245, 250, 251, 398, 409–412, 414–418, 452
- Ligand 2, 6, 11, 14, 15, 17, 18, 20, 21,
23, 51, 52, 54, 55, 61, 63–65, 67–70, 72, 78, 79,
81, 86, 92, 94, 105–107, 109, 111–113, 116, 122,
125, 127, 128, 166, 169, 170, 185, 193, 195, 219,
250, 266, 281, 283, 344, 355, 415, 417, 451
- Linker peptide 18–20, 328
- Lycopersicon esculentum* (tomato) agglutinin
(LEA) 96, 99, 100, 102

M

- Mannose (man) 80, 403, 412
- Mass spectrometry 51–61, 122, 125
- Microcentrifuge 44, 45, 47, 117, 121,
206, 208–210, 224, 225, 334, 347, 424, 453, 456
- Microinjection 328, 330, 333–339
- Monosaccharide 66, 73, 116, 127, 128
- Morpholino (MO) 328–339
- mRNA 328–330, 335, 337–339, 356, 431–433
- MUC (MUC1–MUC17) 37–39, 46, 48, 49, 76, 78
- Mucin 39, 46, 75–89
- Multivalent carbohydrates 80, 81, 84

N

- N-acetylgalactosamine (GalNAc) 76–87
- Neuraminidase 93, 95, 97, 101
- N-glycolylneuraminic acid 79
- N-hydroxysuccinimide (NHS) 65, 69, 95,
97, 118, 121, 124, 271
- Nuclear export 466, 469, 475, 479, 480, 482
- Nuclear magnetic resonance (NMR) 77, 78, 170, 171,
174–176, 183
- Nuclease 329, 331

O

- Oligosaccharide 66, 76, 415
- Oligosaccharide specificity 63
- O-linked 75–77
- Oxidation 17, 18, 52, 53, 60, 61, 350, 412, 415, 416
- Oxidative inactivation 17, 18, 51, 54, 55,
61, 125, 187, 192, 193, 415
- Oxidize 17, 18, 52, 94, 99, 187, 416

P

- Phosphatidylserine (PS) 11, 165, 185–199, 202, 279, 457
- Polymer 19, 80
- Poly-N-acetyllactosamine (polyLacNAc) 19, 20, 71
- Porcine submaxillary mucin (PSM) 77–86, 88
- Precipitate 86, 89, 182, 263, 362, 441, 442, 456
- Pro-angiogenic factors 285, 302
- Prokaryotic 329
- Protein kinase C (PKC) 294
- PS. *See* Phosphatidylserine (PS)
- Pull-down assay 39, 40, 44–48
- Pyridylamination 64, 65

Q

- Quantitative precipitation inhibition 78

R

- Radioactivity 78, 441
- Reducing reagents
- betamercaptoethanol (β ME) 40, 44, 45, 48,
52, 53, 60, 95, 98, 101, 117, 118, 121, 122, 124,
125, 188, 190, 193, 252–254, 259, 260, 426, 428
 - dithiothreitol (DTT) 40, 48, 52, 53, 59,
60, 154–159, 162, 165, 166, 171, 193, 379, 402,
403, 412, 415, 416, 434–439, 441, 444, 445, 447
 - reduced glutathione 52, 53, 171
- Regulatory T cells (Tregs) 5, 16, 251–254, 256–260
- Residue 38, 39, 47, 48, 53, 75–86, 91, 94, 98, 171,
176, 183, 348, 357, 415, 444, 472–474, 478, 481

S

- Saccharide 63, 64, 73, 87, 88, 415,
437–438, 445, 446, 448
- Saturating concentration 93, 195, 276, 278
- SBA. *See* Soybean agglutinin (SBA)
- Sephacrose
- glutathione-conjugated 39, 40, 42–45, 47, 48
 - lactose-conjugated 39, 40, 42–44, 48
- Serine (ser) 16, 76, 78–84, 87, 431, 466, 479
- Sialic acid 20
- Sialylation 20, 73
- Solid phase assay 20, 91–103
- Soluble 21, 52, 61, 76, 147, 157, 180,
270, 281, 282, 293, 295, 299, 306, 375, 460
- Solution 20, 55, 65, 78, 96, 106, 117, 134,
155, 170, 187, 202, 220, 233, 253, 271, 290, 296,
307, 321, 333, 347, 376, 398, 426, 434, 456, 467
- Soybean agglutinin (SBA) 78–88
- Specificity 19, 20, 47, 52, 65, 71,
73, 78, 91, 92, 94, 106, 115–128, 166, 351, 367,
398, 410, 416, 422, 448, 451, 452, 460
- SPR. *See* Surface plasmon resonance (SPR)

Stoichiometry 80, 83, 89
 Subunit 3, 80, 84
 Surface plasmon resonance (SPR) 78, 105–113, 116
 Synthetic polymers 80

T

Tandem repeat 3, 18–20, 47, 76–79, 328
 T cell receptor (TCR) 6, 9, 236, 245, 405
 Terminal glycan 91
 Terminal residue 91
 Tetramer 80, 83–86
 Tetrasaccharide 79, 82
 T helper cells type 2 (Th2) 245, 250, 251
 Thermodynamics 79–83, 116
 Thioether 52, 53
 Threonine (Thr) 76, 78, 79
 TNF. *See* Tumor necrosis factor (TNF)

Transforming growth factor beta (TGF- β) 249
 Translation 38, 125, 249, 328–332, 334, 338, 422, 478
 Tregs. *See* Regulatory T cells (Tregs)
 Tumor metastasis 328
 Tumor necrosis factor (TNF) 9, 14, 21, 142, 148, 185, 218, 297, 300, 454, 457

U

Unfolded protein response 52, 60

V

Vatairea macrocarpa lectin (VML) 79–80, 82–88
 Viral genome 265, 329

Z

Zebrafish 13, 327–339