
Subject Index

- Acetylcholine Receptor (AChR), 69, 71, 88
– muscarinic, 286
– nicotinic, 26, 72
Acetylcholine (see Neurotransmitter)
Actin, 2, 3, 5–8, 12, 52, 54, 59, 71, 90, 133, 135, 179, 182, 186, 219, 220, 222–228, 230–233, 247, 249, 252, 273
Action potential, 129, 130, 134–137, 139, 141
Activated endothelium, 245, 246
Active zone, 4, 49–52, 58, 60–64
– molecular organization, 50, 52, 58, 61
– ultra structure, 51
ADAP, 250
Adaptive Immune Response, 168, 199
Adenosine, 154
Adhesion Molecules, 1, 3, 4–13, 52, 62–64, 70, 74–76, 80, 87, 88, 131–135, 141, 155, 164–167, 176–178, 206, 229, 230, 245, 246, 263, 273, 282
Agrin, 69, 71, 72, 88
AKAP, 8, 12
Alloantigen, 162, 163, 165
Amacrine cell, 115
Anandamide (AEA), 285–287, 289–291, 293
Ankyrin, 2
– AnkyrinG, 132, 133, 138
Anoikis, 158
Antigen Presenting Cell (see also APC), 175–177, 180, 182, 183, 185, 190–193, 199, 200, 202, 203, 205, 207–214, 223, 226, 229–231, 272
AP-1, 289
Apolipoprotein E receptor 2 (ApoER2), 69
Apoptosis, 8, 158, 160, 201, 282, 289
Arachidonic acid, 285, 287
2-Arachidonoylglycerol (2-AG), 286, 287, 290
ARF-GEF, 251, 253
ARNO (Cytohesin-2), 251, 252
Astrocyte, 75–77, 102–104, 134, 138, 264, 269, 282
Autoimmunity, 157, 161
Axoglial junction, 135, 138–141
Axon/axonal, 12, 29, 36, 58, 62, 64, 80, 110, 111, 131–134, 136–144, 283, 284, 287, 288, 290, 291, 293
Axonal Degeneration, 284
B7-H1, 159
B7RP1 (ICOS-L), 158, 159
B Cell Receptor, 202, 203, 213
Basal Ganglia, 107, 288, 290, 291
Basement Membrane, 153, 282
Bassoon, 49, 52, 58–64
Bipolar cell, 74
Blood Cerebrospinal Fluid (CSF) Barrier, 262
Brain Endothelial Cell Lines, 264
Brain Injury, 291
Brevican, 85–88
C2 domain, 49, 53, 55–60
Calcium channel (Ca²⁺ channel), 50, 52–54, 57, 60, 61, 63, 70, 72, 73, 76, 83, 84, 288, 290
– VDCC, 73, 80, 83
Cannabinoid Receptors, 285f.
– CB1, 285, 287–291
– CB2, 285, 289–291
Cannador, 293
Caspr, 88, 132, 134–136, 138–142
CAST/ERC, 49, 52, 58, 59
CC Chemokine Receptor (CCR), 269
– CCR7, 153, 191, 204, 213
– CC2, 269
– CxCR-6, 190
Cdc42, 5–7, 219, 220, 223, 228, 232
Cell adhesion molecule (CAM), 4, 70, 80,

- 87, 133–136, 141
 Ceramide, 289
 Cerebellum, 10, 28, 29, 40, 74, 84, 106, 107,
 111, 112, 117, 133, 288, 290, 291
 Chemokine Ligand (CCL)
 – CCL2, 269
 – CCL19, 153, 213, 264, 269
 – CCL21, 153, 165, 191, 213, 264, 269
 – CXCL12, 153, 227
 – CXCL13, 153
 Chemokine Receptor, 190
 – CC2, 269
 – CXCR-6, 190
 Chemokine/chemokines, 153–155, 165, 168,
 190–192, 204, 206, 213, 227, 245, 246,
 248, 263, 264, 269, 282, 292
 Chondroitin sulfate proteoglycan (CSPG),
 69, 84, 85, 87
 Cluster of Differentiation (CD)
 – CD8, 155–158, 160–163, 165, 166,
 187–190, 192, 201, 211
 – CD28, 158, 179, 232, 233
 – CD31 (PECAM), 164, 248, 273
 – CD34, 153
 – CD40, 158, 187
 – CD80, 158, 159, 179, 186, 233
 – CD86, 158, 159, 186, 233
 – CD95, 282
 – CD98, 248
 – CD99, 164, 273
 – CD154, 158
 Co-Stimulation, 233
 Collagen, 84, 207, 209, 210
 Colloid Osmotic Pressure, 152
 Common Leukocyte Antigen (CLA),
 155
 Complement, 115, 167, 178
 Connexin, 99–107, 110, 113–119
 – Cx36, 99, 101, 104, 106–114, 116–119
 – Cx45, 99, 103–105, 114
 – Cx57, 99, 115
 Connexon, 100
 Contactin, 88, 133–136, 141, 142
 Cross-Presentation, 156, 160
 Cyclooxygenase 2 (COX2), 292
 CYTIP, 252, 253
 Cytohesin, 251–253
 Cytomatrix, 49, 50
 – postsynaptic (see postsynaptic density
 (PSD))
 – presynaptic (see Cytomatrix at the Active
 Zone (CAZ))
 cytomatrix at the active zone (CAZ), 49–52,
 55, 56, 58, 60, 61
 – assembly, 49, 50, 58, 62, 64
 Cytoprotection, 163
 Cytoskeleton, 3, 8, 59, 71, 76, 81, 90, 101,
 133, 135, 219–222, 226, 232, 233, 247,
 251, 252, 273
 Cytotoxic T Lymphocyte (CTL), 157, 162,
 163, 205, 231, 232
 Cytotoxicity, 160, 162, 231, 245, 246
 Death neglect, 201
 Dendrite, 12, 29, 32, 37, 71, 74, 75, 211
 Dendritic Cells (DC), 155–157, 160, 161,
 163, 164, 168, 176, 183, 188, 200, 205,
 248, 259, 267, 274, 284
 Dendritic Spine, 1, 2, 5–13, 75, 78, 79
 Depolarization-Induced Suppression of
 Inhabitation (DSI), 290, 291
 Diacylglycerol (DAG), 49, 50, 55, 56, 177,
 287
 Diapedesis, 166, 263–266, 269, 271–274
 Dictyostelium discoideum, 226
 DNAM, 249
 Dynein (see Motor Protein)
 E-Selectin, 268
 Electrical synapse (see synapse, electric),
 102, 105, 106, 116
 Electron microscopy, 3, 28, 62, 131, 271
 Electroretinogram (ERG), 114
 Electrotonic transmission, 99, 106
 Ena/VASP, 224, 248, 250
 Endocannabinoid System, 281, 285,
 292–294
 Endocytosis, 5, 10, 38–40, 51, 59, 153, 156,
 158, 224, 225
 Endoglycan, 153
 Endomucin, 153
 Endothelial Cell, 151, 153, 154, 156–158,
 160, 162, 163, 165, 166, 252, 264, 266,
 267, 271, 273, 282
 Endothelitis, 157, 162
 Ephrin, 6, 10, 12
 Ephrin receptor (EphB), 6, 12
 Excitatory synapse (see synapse, excitatory)
 Exocytosis, 40, 49–51, 57, 60, 63, 84, 178,
 288

- Experimental autoimmune encephalomyelitis (EAE), 140, 191, 259, 261, 264, 265, 267–273, 283, 284, 292, 293
- Extracellular matrix (ECM), 69–71, 74–77, 79–81, 83–85, 87, 88, 90, 153, 245
- Extracellular Signal-Regulated Kinase (ERK) (see protein kinase), 7, 8, 288, 289
- Extravasation, 192, 245, 246, 261, 272
- Fatty Acid Amide Hydrolase (FAAH), 287, 290
- FERM domain, 247, 249
- Fibronectin (FN), 69, 75, 76, 83, 88, 250
- Filopodia, 62, 220
- Fluorescence resonance electron transfer (FRET), 247
- Focal Adhesion Kinase (FAK, see Protein Kinases), 76, 228, 252, 288
- Formin, 220
- Fyn (see protein kinase), 75, 81, 182, 230, 249, 250
- Fyn binding protein (Fyb), 250
- G-Protein, 226, 266, 269, 288, 289
- G-Protein-Coupled Receptor, 154, 263, 288
- GABA receptors, 25ff
- distribution, 28–30
 - GABAA receptor, 25–39, 81–83
 - GABAB receptor, 82, 83
 - Trafficking, 32ff
- γ -amino butyric acid (GABA) (see Neurotransmitter), 25–40, 69, 81, 82, 85, 288, 290, 291
- GABAergic, 30–32, 69, 81, 82, 84, 89, 109, 111, 112, 118, 290, 291
- γ -frequency (gamma-frequency), 110, 111, 116
- Gap junction, 72, 99–102, 104–107, 109, 110, 113, 116–119
- Gephyrin, 29–32, 38, 40
- Glia, 84, 85, 102, 129, 131, 137, 142, 143, 262, 267
- Gliomedin, 133, 134, 138
- Glutamate (see Neurotransmitter)
- glutamatergic, 60, 84, 288, 291
- Glutamate receptor, 1ff, 35, 69, 73, 84, 119, 286, 287
- AMPA receptor, 5–13, 35, 38, 73, 76, 79, 80, 84, 89
 - GluR1, 8, 9, 11, 69, 79
 - GluR2, 5, 9–12, 84
 - metabotropic glutamate receptor, 10, 286
 - NMDA receptor/NR2A/NR2B, 3, 5–11, 70, 73, 75, 76, 78–80, 83, 88, 89, 119
- Glutamatergic, 60, 84, 288, 291
- Glycine receptor, 26, 31, 38
- Golgi apparatus, 33, 35, 37, 59, 63, 64, 178, 229, 231, 232, 251
- GRASP, 10, 252
- GRIP, 5, 9–13, 35, 37
- GRP1, 251, 252
- GTPase, 2, 7, 52, 56, 223, 225, 247, 248
- hDlg (see human Discs large)
- Heart, 155
- Helper T Cells, 190
- Hemichannel, 100
- Hemisynapse, 175, 177, 180, 182, 183
- Heparan sulphate proteoglycan (HSPG), 71, 80
- Heparin-binding growth-associated molecule (HB-GAM), 69, 80, 81
- Hepatocytes, 153
- High Endothelial Venules (HEV), 153, 175, 187–189, 203, 245, 246
- Hippocampal/Hippocampus, 10, 28, 29, 31, 32, 35, 36, 38, 60, 62, 69, 72, 74–76, 78–88, 99, 107, 109–111, 116–118, 287, 288, 290–292
- CA1, 28, 57, 74, 76, 78, 80–86, 88, 104, 109–111, 287, 288, 290
 - CA3, 84, 109, 111
- HNK-1, 81, 82
- Horizontal cell, 74, 115
- Human discs large (hDlg) (see also SAP97), 249
- Hyaluronic acid (HA), 69, 85
- Immune privilege, 259, 260, 274, 281, 283
- Immunological Synapse (IS), 175–193, 207, 208, 219, 224, 229, 230, 231, 271
- Immunoreceptor Tyrosine-Based Activation Motifs (ITAM), 177, 230
- Inflammation, 158, 161, 164, 165, 167, 168, 176, 183, 186, 245, 246, 267, 269–271, 274, 283, 284, 292
- Inhibitory synapse (see synapse, inhibitory)
- Inner ear hair cell, 61
- Innexin, 116, 117

- Inside-out, 77, 246, 250, 251, 263
 Integrin, 69, 70, 73, 75–77, 167, 177, 182,
 219, 230, 232, 245–251, 253, 263, 266,
 268, 270–272
 – $\alpha 4\beta 7$, 263, 270, 271
 Intercellular Adhesion Molecule-1 (see also
 ICAM), 178, 179, 182, 208, 245–248, 250,
 264, 270, 271, 273
 Interferon (IFN), 154, 156, 159, 160, 162,
 164–166, 168, 289
 Interleukin (IL)
 – IL-1 β , 289
 – IL-12, 158, 159, 164, 249, 282, 289
 – IL-2, 160, 162, 164, 212, 224, 231, 248, 249
 – IL-4, 160, 167, 178, 185, 231, 289
 Interneuron, 28, 79, 81, 99, 108–113, 116,
 118, 288, 290, 291
 Internode, 131, 135
 Intravital (Video) Microscopy (IVM), 166,
 185–187, 189, 190, 265, 266, 268, 272
 Ion channel, 25, 26, 31, 69, 131, 134, 135,
 141, 249

 JAB-1, 253
 JAK, 252, 282
 Junctional Adhesion Molecules (JAMs),
 152, 273
 Juxtaparanode, 130–132, 135, 136, 138, 142,
 143

 Kidney, 155, 156
 Kinase (see protein kinase)
 Kinesin (see also motor protein), 12, 13, 54,
 249

 L-Selectin, 188, 270
 Lamellipodia, 182, 220, 225–227, 249
 Laminin, 72–75, 84, 89
 Leading Edge, 138, 180, 193, 206, 221, 227,
 248
 Lecticans, 85, 87
 Leucocyte Function-Associated Antigen
 (LFA)
 – LFA-1, 154, 165, 178, 182, 183, 192,
 229–232, 245, 246, 248–250, 253, 264,
 266, 270, 271, 273
 – LFA-3, 158
 Linker for activation of T cells (LAT), 177,
 230–232
 Lipopolysaccharide (LPS), 188, 245, 246,
 248, 260, 268
 Liprin, 9, 10, 13, 52, 58, 64
 Live Cell Imaging, 208, 212
 Liver, 152, 153, 155–157, 160, 161, 163, 167,
 168, 177, 190, 191
 Liver Sinusoidal Endothelial Cells (LSEC),
 152, 153, 157, 163, 168
 Long-term depression (LTD), 8–10, 83, 87,
 88, 291
 Long-term potentiation (LTP), 8–11, 57, 69,
 70, 74–76, 78–83, 85–88, 90, 290
 Lung Microvascular Endothelial Cells, 152,
 163
 Lymph Node, 153, 156, 176, 184–188,
 190–192, 200, 202–205, 208, 210–213,
 248
 Lymphatic Endothelial Cells, 153, 168

 Macrophage, 156, 157, 200, 201, 203, 208,
 224, 228, 274, 282–284, 289, 291
 Macropinocytosis, 156
 Major Histocompatibility Complex (MHC),
 155–157, 159, 160, 163–165, 176–186,
 192, 200–203, 205, 208, 210, 229, 230,
 259, 282
 MAP kinase (see protein kinase/ERK)
 Marinol, 293
 Matrix metalloprotease (MMP), 70, 88, 90
 Membrane-associated guanylate kinase
 (MAGUK) (see also SAP97, PSD-95,
 hDlg), 249
 Memory Cells, 206
 Microenvironment, 87, 90, 151, 153, 158,
 259, 264, 282
 Microglial Cells, 177, 192, 267, 282–285,
 290–292
 Microtubule (see tubulin), 12, 31, 34, 178,
 207, 219, 221–223, 226, 229–231, 249
 Microtubule Organizing Centre (MTOC),
 207, 221, 222, 227–232
 Migration
 – of neurons, 70, 76, 77
 – of leucocytes, 164–167, 175, 176, 183–187,
 189–191, 200, 202, 203, 206, 210–213,
 219, 221, 223, 226, 227, 229, 248, 250,
 253, 259, 261, 263–270, 273, 274
 Minocycline, 292
 Mitogen-activated protein kinase (MAPK)
 (see protein kinase), 8, 72, 178
 Monoacylglycerol Lipase (MGL), 287

- Monocyte Chemoattractant Proteins (MCP-1), 228
- Motor protein, 1, 12, 61, 249
- Dynein, 232
 - Kinesin/KIF, 12, 13, 249
- Multiple Sclerosis, 140, 259, 260, 265, 283
- Munc13 (UNC13), 50, 52, 53, 55–57, 60, 61, 63, 64, 253
- Muscle-specific receptor tyrosine kinase (MUSK), 70
- Myelin, 102, 103, 129–131, 134, 137, 139–141, 143, 175, 191, 265, 269, 283, 284
- N-Acyltransferase (NAT), 285
- N-WASP, 220, 223–225
- Narp (pentraxin), 69, 70, 79
- Natalizumab, 270, 271
- Natural Killer (T) Cells, 81, 177
- Naïve CD4 T Cells, 157, 160, 164
- NCAM, 64, 73, 80, 82, 87
- Neocortex, 75, 106, 108, 288, 290
- Neurodegenerative Diseases, 283, 292, 293
- Neurofascin, 132, 133, 135
- Neurogenesis, 101, 285
- Neuroinflammation, 284, 293
- Neurological disorder/disease, 281
- Neuromuscular junction (NMJ), 51, 69–72, 74, 77, 88
- Neuron-Glia, 129, 131–133, 136–144
- synapse (glia-neuron), 84
- Neuroprotection, 281, 283, 293, 294
- Neurotransmitter, 25, 31, 49–51, 56–58, 60, 118, 282, 288, 290, 291
- acetylcholine, 26, 69–71, 286, 291
 - γ -aminobutyric acid (see also GABA), 25–40, 290, 291
 - dopamine, 10, 115, 118, 119, 286, 290
 - glutamate, 5, 75, 84, 286–288, 290, 291
 - release, 50, 51, 56–58, 60–62, 73, 82, 87, 99, 118, 253, 282, 288, 291
- NF-AT, 289
- NG2, 84, 85
- Nitric oxide (NO), 7, 115, 116, 154, 289, 292
- Nitric Oxide Synthase (NOS), 7, 79, 85, 90, 99, 101, 102, 104–110, 113, 116–118, 129, 131, 223, 260, 264, 281–288, 291–293
- Node of Ranvier, 129–131, 133, 134, 136, 138, 140, 143, 144
- Non-Steroidal Anti-Inflammatory Drugs (NSAID), 292
- Oligodendrocyte, 84, 102, 103, 129, 130, 134, 138, 139, 140, 144, 269
- Oral Tolerance, 161, 185
- p38 Mitogen-Activated Protein Kinase (see protein kinase)
- Pain, 291, 293, 294
- Palmitoylation, 1, 10, 11, 37
- Pannexin, 99, 100, 106, 116–118
- Paranode, 131, 134–136, 140–144
- PDZ domain, 3, 4, 10, 11, 30, 35, 57, 58, 60, 64, 143, 252
- PDZ-binding motif, 7, 249
- PECAM (see CD31), 164, 248, 273
- Peyer's Patches, 153
- Phage Display Analysis, 151
- Phospholipase D, 285
- Photoreceptor, 49, 59–61, 74
- Piccolo, 49, 52, 58–64
- Piccolo-Bassoon transport vesicle (PTV), 49, 62–64
- Pioglitazone, 292
- Plasticity (see synaptic plasticity), 1, 5, 7–9, 13, 40, 56, 57, 59, 69, 71, 74–76, 78, 80, 81, 83–85, 87, 88, 90, 294
- PMA, 248, 249, 252
- Podocalyxin, 153
- Podosomes, 223
- Polysialic acid (see NCAM), 64, 80, 82, 87
- Postsynaptic density (PSD), 1–7, 9–13, 49–51, 56, 136
- Potassium (K^+) channels, 5, 11, 69, 80, 82, 83, 132, 135, 136, 142–144, 249, 288, 289
- Kv, 136, 143
- Presynaptic bouton terminal, 30, 32, 51, 56, 61, 71, 72, 82, 90, 287, 288, 290, 291
- Profilin, 248
- ProSAP (also Shank), 2–4, 6, 7, 12
- Prostacyclins, 154
- Proteasome, 11, 35, 36, 157
- Protein kinase
- CaMKII, 7, 8, 11, 72
 - ERK, 7, 8, 288, 289
 - FAK, 76, 228, 252, 288
 - Fyn, 81, 182, 230, 249, 250
 - p38, 8, 288
 - PKA/protein kinase A/cAMP-dependent

- protein kinase, 8, 36, 57, 60, 115, 285, 288
 – Src, 8, 72, 76, 78, 182, 230, 250
 – ZAP-70, 177, 179, 181, 182, 230–232
 PSD-95, 2–12, 30, 32, 37, 73
 Purkinje cell, 74

 Rab3, 49, 52, 56, 57, 59, 63
 Rac, 2, 220, 225
 Rap1, 219, 247–249
 Ras, 2, 7, 8, 10, 177, 219
 Receptor protein tyrosine phosphatase (RPTP), 9, 70, 87, 88
 Reelin, 69, 73, 77–79
 Regulatory T Cells, 161, 183, 206
 Reticulo-Endothelial Cell System, 156
 Retina, 4, 31, 59, 61, 74, 77, 99, 101, 107, 110, 113–115, 119, 139
 Rho, 5, 7, 88, 133, 219, 220, 223, 225, 226, 286
 RIAM, 248, 249
 Ribbon synapse, 49, 59–61
 RIM, 50, 52, 56–60, 63, 64

 SAP97, 2–4, 8, 9, 11, 12 (see also hDlg)
 Sativex, 293
 Scaffold, 8, 25, 31, 32, 37, 55, 58, 61, 63, 71, 73, 79, 83, 133, 138, 252
 Schwann cell, 132–134, 137, 138, 144
 SDF-1, 227
 Selectins, 154, 206, 268
 Selection, 201, 202
 Serial Encounter Model, 199, 210, 212
 SH2 Domain-Containing,
 – Leucocyte Phosphoprotein of 76 kDa (SLP-76), 181, 230–232, 250
 Shank, 2–4, 6, 7, 12
 Shiverer Mouse Mutant, 140
 Sialomucins, 153
 Signaling proteins, 1, 3, 11, 245, 246
 Single Encounter Model, 199, 208
 siRNA/RNAi, 5, 7–9, 12, 32, 85, 221, 225, 234, 247, 250
 SKAP-55, 250
 Skin, 117, 155, 167, 176, 203, 205
 Sodium (Na⁺) channels, 129–133, 137–144
 Spectrin, 133, 135
 Spectrin/Fodrin, 7, 53, 133, 135
 Sphingomyelin (SM), 50, 289
 Spine (dendritic spine), 5–8, 13, 78
 Stargazin, 5, 9, 11, 73, 79
 Superantigen, 222

 Supra-Molecular Activation Clusters (SMAC), 178, 182, 230–232
 Synapse (see also immunological synapse),
 1ff, 25ff, 49ff, 69ff, 99ff, 129, 175ff, 230, 287–291
 – chemical, 49, 50, 72, 86, 90, 99, 101, 102, 105, 106, 116, 118, 119
 – electrical, 99ff, 101, 102, 105–108, 110–113, 115, 116, 118, 119
 – excitatory, 1ff (1–13), 30–32, 57, 69, 75–79, 81, 291
 – formation/developmental/assembly, 3, 4, 7, 11, 12, 30, 33, 49, 50, 56, 58, 62–64, 69, 72, 75, 77, 84, 85, 89, 101
 – history, 1
 – immunological (see immunological synapse)
 – inhibitory, 4, 5, 25ff, (25–40), 57, 81, 82, 86, 291
 – mixed, 119
 Synaptic plasticity, 5, 7–9, 40, 59, 69, 71, 73–78, 80, 81, 83, 87–89, 294
 Synaptic transmission, 1, 9, 10, 25ff, 56, 57, 61, 69, 72, 74, 81–83, 85–87, 90, 99, 106, 116, 118, 119, 219, 271, 291
 Synaptic vesicle (SV), 30, 49–52, 55–57, 59, 61–63, 70, 74, 82
 Synaptogenesis (see also synapse formation), 30, 35, 49, 50, 58, 64, 69–72, 74–77, 79, 80
 Syndecan, 7, 80, 81
 SynGAP, 7, 8

 T Cell,
 – encephalitogenic, 265–271, 273
 – Receptor, 165, 193, 200, 209, 224, 229, 245, 246
 – Tolerance, 160, 161, 168, 176, 177, 183–187, 189, 190, 193, 211
 Talin, 178, 247
 Tenascin, 69, 70, 81, 83, 134
 – Tenascin-C (TN-C), 83, 84
 – Tenascin-R (TN-R), 81–83
 Tetrahydrocannabinol (THC), 285, 289, 293, 294
 Th1 cells, 160, 164, 249
 Thrombospondin (TSP), 70, 77
 Thymic Organ Cultures, 210
 Tight Junction (TJ), 101, 152, 260, 272
 Time-Lapse Microscopy, 199

- Two Photon Laser Scanning Microscopy (TPLSM), 184, 185, 187
- Trafficking,
– GABA-receptor, 25, 32–37
– glutamate receptor, 1, 6–13
– leukocyte, 153, 245, 246, 261, 264, 265, 269–271, 273
– membrane, 49–54, 59
– proteins, 3
- Transcytosis, 152, 153
- Transendothelial migration, 154, 166, 252
- Transendothelial Migration (Transmigration), 154, 164–168, 252, 263
- Transport,
– axonal, 13, 49, 54
– neurotransmitter receptors, 4, 12, 34, 39
– vesicular, 37, 49, 62–64
- Transporter,
– dopamin, 10
- Tumor necrosis factor (TNF), 154, 156, 159, 160, 167, 176, 260, 268, 282, 289
- Ubiquitin/ubiquitinylation, 11, 12, 35, 36
- Uropod, 180, 183, 191, 206, 227, 228
- VASP, 224
- Vav1, 230–233
- VCAM-1 (Vascular Cell Adhesion Molecule-1), 263, 264, 266, 270, 272
- Very low density lipoprotein (VLDL), 70, 78
- Vessels, 151, 153, 154, 161, 162, 202, 203, 206, 259, 262, 265–270
- WASP, 220, 221, 223–225, 227, 228, 230, 232
- WASP Interacting Protein WIP (WASPIP), 223, 224
- WAVE, 225–228
- Wiskott Aldrich Syndrome (WAS), 223, 228
- ZAP-70 (see Protein Kinase), 177, 179, 181, 182, 230–232
- ZO-1, 49, 56, 152
- ZO-2, 152